

Review of ONT Sequencing Runs

Dr. Louise Judd
Centre for Pathogen Genomics- Innovation Hub
Lead Laboratory Scientist and Genomics Trainer


How do we do genomic sequencing?

1. Clinical sample




- High viral load
- Skin swabs
- Swab storage
- DNA extraction
- Negative extraction controls

2. Initial test to see if there's a bug



- Q-PCR
- High viral load
- Ct <29
- DNA storage

3. Make more copies of RNA or DNA



- Tiled amplicon PCR
- Artic 2500 bp scheme


4. Sequencing



- ONT platform
- Rapid libraries
- Sup or HAC basecalling
- 50K reads at N50 ~1200 bp

How do we do genomic sequencing (cont.)?

5. Put the pieces together




Generate
consensus
genome

6. Line them up

```
AGTCGATTAGCTTAGCTTGTAGCGCTATATTAT
AGTCGATTAGCTTAGATTGTAGC-CTATATTAT
AGTCGATTAGCTTAGATTGTAGC-CTATATTAT
AGTCGATTAGCTTAGATTGTAGC-CTATATTAT
AGTCTTATTGCTGAGCTTGTAGC-CTAT--AT
AGTCGATTAGCTTAGCTTGTAGCGCTTATTAT
AGTCGATTAGCTTAGCTTGTAGCGCTATATTAT
AGTCGATTAGCTTAGCTTGTAGCGCTATATTAT
AGCTGTTTAGCTTAGCTTGTAGCGCTATATTAT
```

Align all the sequences
to reference genome to
compare

7. Map the similarities
and differences



Visualize relationships
between sequences

8. Put it all together



Mpox genomics recap


Source Details	Barcode	PCR DNA pool 1 ng/ul	PCR DNA pool 2 ng/ul	final library ng/ul
sample1	barcode41	7	16	
sample2	barcode42	41	46	
sample3	barcode43	5	11	
sample4	barcode44	23	26	
sample5	barcode45	11	11	
sample6	barcode46	30	45	
sample7	barcode47	5	5	
NTC	barcode48	6	12	38

Where do I find my sequencing data?

1. Default location for run is “data”

2. Run name will be as entered when the run was being started

3. Will by default say no_sample_id unless sample ids used at set up



New folder each time run is started.

Folder name

- Date stamp
- Time stamp
- Sequencing device ID (MN)
- Flow cell ID
- Random unique number

How is my sequencing data structured?

Name	Date modified	Type	Size
fastq_fail	28/05/2025 ...	File folder	
fastq_pass	28/05/2025 ...	File folder	
other_reports	27/05/2025 ...	File folder	
pod5	28/05/2025 ...	File folder	
barcode_alignment_FBB...	28/05/2025 ...	TSV File	13 KB
final_summary_FBB3034...	28/05/2025 ...	Text Docu...	1 KB
pore_activity_FBB30344...	28/05/2025 ...	Excel.CSV	372 KB
report_FBB30344_2025...	28/05/2025 ...	Microsoft ...	1,324...
report_FBB30344_2025...	28/05/2025 ...	JSON Sou...	3,928...
report_FBB30344_2025...	28/05/2025 ...	Markdow...	453 KB
sample_sheet_FBB30344...	28/05/2025 ...	Excel.CSV	1 KB
sequencing_summary_F...	28/05/2025 ...	Text Docu...	522,1...
throughput_FBB30344_...	28/05/2025 ...	Excel.CSV	80 KB

fastq_fail

- reads with Q score below cut off (<10)
- Discard

fastq_pass

- reads with Q score above cut off (<10)
- One directory per sample or barcode

pod5

- raw data files
- very large files
- Only needed if you wish to rebsecall
- Keep until basecalling successful

reports

- Move all files to other_reports
- Save whole directory

Sequencing run report

- html file that can be opened with any web browser
- Will only be made if the run finishes successfully

GridION GXB04127 Final report

17 Jun 25, 15:21 UTC+10:00 — 18 Jun 25, 09:38 UTC+10:00 · O2025-033_MPOX_C · O2025-033 · X2

Protocol run ID: aab53bbb-8364-459c-9c81-e05cf536db81

[Run summary](#) | [Run configuration](#) | [Sequence output](#) | [Run health](#) | [Run log](#)



Sequencing instrument

Date and time stamp of run

Run name

Unique run ID

Sequencing run report- Run summary

^ Run summary

DATA OUTPUT

Estimated bases

3.25 Gb

Reads generated

3.39 M

Estimated N50

1.2 kb

BASECALLING

Reads called

100%

Reads called (min Q score: 10)

2.81 M

687.53 k

Pass

Fail

Bases called (min Q score: 10)

1.86 Gb

313.57 Mb

Pass

Fail

RUN DURATION

Run time

18 hrs 12 mins / 24 hrs 0 mins

Elapsed time



Run limit



Run status

STOPPED · By user

Check basecalling is at 100%

N50: 1.0-1.5 kb

Proportion reads failing Q Score <10

Sequencing run report- Run configuration

Run configuration

Run Setup

Flow cell type	FLO-MIN114
Flow cell type alias	FLO-MIN114
Flow cell ID	FBA40009
Kit type	SQK-RBK114-96

Run Settings

Run limit	24 hrs
Pore scan freq.	1.5 hrs
Reserved pores	On
Basecalling	Super-accurate basecalling v4.3.0, 400 bps
Modified basecalling	Off
Trim barcodes	On
Mid-read barcode filtering	Off
Min Q score	10

Data Output Settings

FAST5 output	Off
FASTQ data output	One file per hour
POD5 data output	One file per hour, or 500000000 bases per batch
BAM file output	Off
Bulk file output	Off
Data location	/data/./O2025- 033_MPOX_C/O2025- 033/20250617_1521_X2_F BA40009_aab53bbb

Software Versions

MinKNOW	24.11.8
Bream	8.2.5
Configuration	6.2.12
Dorado	7.6.7
MinKNOW Core	6.2.6

Sequencing run report- Read length distribution

Sequence output

READ LENGTHS · OUTLIERS REMOVED


The read length graph shows the total number of bases vs the read length. The longest 1% of strands are classified as outliers, and excluded to allow focus on the main body of data.

N50*

1.2 kb

% Basecalled

100 %



OUTLIERS

The longest 1% of strands are classified as outliers, and aggregated into groups to show their relative amounts.

Read length (kb)	Bases (Mb)
0 - 131.072	17.82
131.072 - 262.144	0.44
262.144 - 389.12	0.39

*N50 calculated from basecalled read length histogram.




Sequencing run report- Base output


^ CUMULATIVE OUTPUT

The cumulative output shows the total amount of bases or reads sequenced over time by your device.

Bases

Legend

 Estimated	 Passed	 Failed
Predicted total number of bases, prior to basecalling	Bases equal to or above the quality score threshold.	Bases below the quality score threshold.



Sequencing run report- Read output

Reads

Legend

— Total

Total number of reads, including passed, failed and skipped.

— Passed


Reads equal to or above the quality score threshold.

— Failed

Reads below the quality score threshold.

— Skipped

Reads that will not be basecalled. Post run basecalling is possible.



Sequencing run report- Barcode distribution table

BARCODES

Detected barcodes Bases graph Reads graph

Detected barcodes

The total number of bases and reads for each barcode detected are displayed in table below. Reads/bases must have a quality score above 9 to pass.

Search barcodes

Q

Unclassified data

133.22 k (4.7%)

Passed Reads

103.62 Mb (5.6%)

Passed Bases

Export CSV



Barcode	Total bases (Mb)	Passed bases (%)	Total reads (k)	Passed reads (%)
barcode01	58.9	93	102.53	88.6
barcode02	83.56	93.4	131.52	90
barcode03	1.31	44.3	15.98	45.1
barcode04	52.79	90.1	104.59	85.6
barcode05	12.16	88.9	41.92	80.8
barcode06	68.86	92.1	126.73	87.1
barcode07	1.72	56	25.7	53.5
barcode08	1.22	51.9	14.81	48.4
barcode09	63.03	91.3	75.84	90
barcode10	76	91	84.27	90.1

Barcode distribution csv file- open with Excel

AutoSave

O2025-033_MPOX_C_FBA40009_barcode (1)

Home

Insert

Draw

Page Layout

Formulas

Data

Review

View

Automate

Paste

Aptos Narrow (Bod... 12 A^ A^

B I U

General

Conditional Form

Format as Table

Cell Styles

G9

	A	B	C	D	E	F	G	H	I	J	K	L
1	Barcode	Total bases (Mb)	Passed bases (%)	Total reads (k)	Passed reads (%)							
2	barcode01	58.9	93	102.53	88.6							
3	barcode02	83.56	93.4	131.52	90							
4	barcode03	1.31	44.3	15.98	45.1							
5	barcode04	52.79	90.1	104.59	85.6							
6	barcode05	12.16	88.9	41.92	80.8							
7	barcode06	68.86	92.1	126.73	87.1							
8	barcode07	1.72	56	25.7	53.5							
9	barcode08	1.22	51.9	14.81	48.4							
10	barcode09	63.03	91.3	75.84	90							
11	barcode10	76	91	84.27	90.1							
12	barcode11	0.68	31.5	6.31	53.3							
13	barcode12	59.95	92	71.95	91.1							
14	barcode13	16.64	93.1	43.4	87.1							
15	barcode14	102.26	91.7	114.41	90.9							
16	barcode15	0.72	51.4	9.93	48.8							
17	barcode16	0.4	42.6	5.42	50.9							
18	barcode17	66.31	89.9	125.44	84.1							
19	barcode18	88.32	92.6	150.11	88.5							
20	barcode19	1.96	42.6	19.36	52.5							
21	barcode20	37.57	91.1	82.89	84.6							
22	barcode21	12.63	86.2	55.71	72							
23	barcode22	79.93	91.4	167.56	84.6							
24	barcode23	1.06	47.9	16.41	51							
25	barcode24	2.86	59.8	29.25	53.1							
26	barcode25	93.66	90.3	130.96	86.4							
27	barcode26	158.78	91.7	190.53	89.9							
28	barcode27	0.45	50.3	7.34	46.7							
29	barcode28	111.94	91.2	147.06	89.1							

Sequencing run report- Barcode distribution graph

^ BARCODES

Detected barcodes **Bases graph** Reads graph

Bases graph

The total number of bases for each barcode detected and displayed in the graph below. Bases must have a quality score above 9 to pass.

Unclassified data

103.62 Mb


Passed Bases

5.6%

Percentage

Barcodes ascending ▾

Expand ⬆




Sequencing run report- Q Score distribution graph

^ QUALITY SCORE

The quality score is calculated as basecalling is performed on your device. Reads that fall below the minimum value of 9 will be classified as failed reads. You can alter the accepted minimum quality score in MinKNOW.

Q Score histogram

Passed simplex bases Failed simplex bases




Filters

Min Q score 10

Q Score over time

Legend

- Mode
The most frequent quality score of reads in the run.
- Spread
The spread of quality scores, found by calculating full width half maximum.
- - - Min. quality score
Minimum quality score to be accepted as a passed read.



Sequencing run report- Pore Activity

PORE ACTIVITY


The Pore activity graph shows the performance of your sample as it is being sequenced during a run.

☒ Show grouped

Legend

- Sequencing
Pore currently sequencing
- Pore available
Pore available for sequencing
- Unavailable
Pore currently unavailable for sequencing
- Inactive
Pore no longer suitable for further sequencing
- Unclassified
Pore status unknown

Aim: Pore availability to stay constant for first 24 hours of run



Sequencing run report- Pore Scan

PORE SCAN

A Pore scan is performed at configurable time intervals to determine the current status of pores within channels on a Flow Cell. For this run a Pore scan is performed every 1.5 hrs.

Legend

- Pore available

Pore in channel available for sequencing
- Reserved pore


Pore in reserve, will return to available when required
- Unavailable

Pore inhibited from sequencing
- Saturated





Possible contamination in the sample
- Zero

No current is passing through this pore, possibly due to bubbles on the membrane
- Inactive

Pore no longer suitable for further sequencing



Unavailable pores:
Can become
available again
after nuclease
wash




Washing and reusing flow cells

PORE SCAN

A Pore scan is performed at configurable time intervals to determine the current status of pores within channels on a Flow Cell. For this run a Pore scan is performed every 1.5 hrs.

Legend

- Pore available
Pore in channel available for sequencing
- Reserved pore
Pore in reserve, will return to available when required
- Unavailable
Pore inhibited from sequencing
- Saturated
Possible contamination in the sample
- Zero
No current is passing through this pore, possibly due to bubbles on the membrane
- Inactive
Pore no longer suitable for further sequencing



Flow Cell Wash Kit (EXP-WSH004 or EXP-WSH004-XL)



Version: WFC_9120_v1_revQ_08Dec2020
Last update: 4/25/2024

Kit batch number Flow cell number DNA Samples

Checklist: Flushing a MinION/GridION Flow Cell


Materials	Consumables	Equipment
<input type="checkbox"/> Flow Cell Wash Kit (EXP-WSH004) or Flow Cell Wash Kit XL (EXP-WSH004-XL)		<input type="checkbox"/> P1000 pipette and tips <input type="checkbox"/> P20 pipette and tips <input type="checkbox"/> Ice bucket with ice

TRANSLOCATION SPEED

The translocation speed is the rate at which DNA/RNA travels through pores as it is being sequenced.

Legend

- Median
- 75% quartile
- 25% quartile
- Accepted range




TEMPERATURE

The temperature of the Flow Cell over the run time.

Legend

- Measured
- Target



Mpox genomics recap

Source Details	Barcode	PCR DNA pool 1 ng/ul	PCR DNA pool 2 ng/ul	final library ng/ul
sample1	barcode41	7	16	
sample2	barcode42	41	46	
sample3	barcode43	5	11	
sample4	barcode44	23	26	
sample5	barcode45	11	11	
sample6	barcode46	30	45	
sample7	barcode47	5	5	
NTC	barcode48	6	12	38

Seqkit Stats to generate read stats

<https://bioinf.shenwei.me/seqkit/>

simple statistics of FASTA/Q files

Columns:

- | | |
|-------------|---|
| 1. file | input file, "-" for STDIN |
| 2. format | FASTA or FASTQ |
| 3. type | DNA, RNA, Protein or Unlimit |
| 4. num_seqs | number of sequences |
| 5. sum_len | number of bases or residues , with gaps or spaces counted |
| 6. min_len | minimal sequence length , with gaps or spaces counted |
| 7. avg_len | average sequence length , with gaps or spaces counted |
| 8. max_len | miximal sequence length , with gaps or spaces counted |
| 9. Q1 | first quartile of sequence length , with gaps or spaces counted |
| 10. Q2 | median of sequence length , with gaps or spaces counted |
| 11. Q3 | third quartile of sequence length , with gaps or spaces counted |
| 12. sum_gap | number of gaps |
| 13. N50 | N50. https://en.wikipedia.org/wiki/N50,_L50,_and_related_statistics#N50 |
| 14. N50_num | N50_num or L50. https://en.wikipedia.org/wiki/N50,_L50,_and_related_statistics#L50 |
| 15. Q20(%) | percentage of bases with the quality score greater than 20 |
| 16. Q30(%) | percentage of bases with the quality score greater than 30 |
| 17. AvgQual | average quality |
| 18. GC(%) | percentage of GC content |
| 19. sum_n | number of ambitious letters (N, n, X, x) |

Source Details	PCR DNA pool 1 ng/ul	PCR DNA pool 2 ng/ul	final library ng/ul	QPCR Ct values	num_seqs	sum_len	N50	Q20(%)	Q30(%)	depth
sample1	7	16		25.6	40,871	23,601,986	1047	83	74.11	118.01
sample2	41	46		21.9	42,570	30,341,778	1233	83.04	74.39	151.71
sample3	5	11		Not detected	6,365	2,029,810	830	76.62	62.16	10.15
sample4	23	26		21	69,823	47,729,766	1245	84.77	76.02	238.65
sample5	11	11		29.1	47,565	15,985,503	529	83	70.01	79.93
sample6	30	45		19	71,881	52,654,060	1365	81.48	72.95	263.27
sample7	5	5		Not detected	3,918	518,858	425	55.61	39.69	2.59
NTC	6	12	38	N/A	8,848	3,008,851	808	74.93	62.26	15.04

Source Details	PCR DNA pool 1 ng/ul	PCR DNA pool 2 ng/ul	final library ng/ul	QPCR Ct values	Kraken Hit 1	% reads	Kraken Hit 2	% reads	Kraken Hit 3	% reads
sample1	7	16		25.6	Monkeypox virus	41.8	Photobacterium leiognathi	0.6	Vaccinia virus	0.1
sample2	41	46		21.9	Monkeypox virus	64.4	Homo sapiens	1.6	Photobacterium leiognathi	0.6
sample3	5	11		Not detected	Photobacterium leiognathi	1.4	Homo sapiens	0.1	Streptomyces sp. NBC_01420	0.0
sample4	23	26		21	Monkeypox virus	79.7	Homo sapiens	2.8	Vaccinia virus	0.3
sample5	11	11		29.1	Homo sapiens	46.4	Monkeypox virus	6.7	Streptomyces sp. NBC_01353	3.1
sample6	30	45		19	Monkeypox virus	80.2	Homo sapiens	4.1	Photobacterium leiognathi	0.5
sample7	5	5		Not detected	Photobacterium leiognathi	0.4	Monkeypox virus	0.0	Escherichia coli	0.0
NTC	6	12	38	N/A	Photobacterium leiognathi	0.6	Homo sapiens	0.1	Streptomyces sp. NBC_01420	0.0

Mpox primers to take home

- ARTIC INRB pool1 and pool2
- 100 uM stock of pooled primers
- 30 ul volume
- PLEASE store them at -20°C freezer as soon as possible
 - In the freezer in your hotel room
 - In your freezer at home
 - In your freezer at the lab

		Concentration ng/ul					
	Mpox QPCR Ct	pool1 from training	pool2 from training	pool1 - 20oC	pool2 - 20oC	pool1 RT O/N then - 20oC	pool2 RT O/N then - 20oC
sample1	25.6	7	16	12	25	4	18
sample2	21.9	41	46	44	50	23	44
sample3	Not detected	5	11	3	4	2	3
sample4	21	23	26	36	38	29	35
sample5	29.1	11	11	9	9	7	8
sample6	19	30	45	43	49	42	43
sample7	Not detected	5	5	3	4	4	2
NTC	N/A	6	12	3	4	4	2